

DeltaLog10 release notes

Revision Number	Date	Notes
1.1.0.6	22/11/2023	Cumulative update that includes updated components and corrections in the "Discomfort Analysis" module. Added Visual C++ redistributable 2013 installation
1.1.0.4	30/05/2023	Added printing of the report for "Discomfort" in moderate environments according to the ISO7730 standard.
1.1.0.2	25/10/2022	Calculation correction in "Cold environments" module.
1.1.0.1	30/09/2022	Improved the connection with the instruments.
1.1.0.0	12/07/2022	General review.
1.0.1.7	14/02/2022	Added CP210x_VCP USB drivers necessary for HD32.1 manufactured from 02/2022.
	18/02/2020	Added compatibility for HD32.3TC.
1.0.1.6	11/09/2019	WBGT lim visualization according to new revision of ISO7243 standard.
1.0.1.5	22/05/2019	Changed default installation folder.
1.0.1.4	30/01/2019	HD21ABE17 monitor bug fixing.
1.0.1.3	20/12/2018	Bugs fixing.
1.0.1.2	13/11/2018	Regenerated setup file.
1.0.1.1	09/11/2018	Added new functionalities to comply with ISO7243:2017 standard. Corrected a calculation inaccuracy of the ISO7933 PHS index.
0.1.6.7	17/10/2018	Added setting of CAV parameter for WBGT correction according to ISO 7243:2017 standard.