

DeltaLog10 release notes

Revision Number	Date	Notes
1.0.1.7	14/02/2022	Added CP210x_VCP USB drivers necessary for HD32.1 manufactured from 02/2022.
	18/02/2020	Added compatibility for HD32.3TC.
1.0.1.6	11/09/2019	WBGT lim visualization according to new revision of ISO7243 standard.
1.0.1.5	22/05/2019	Changed default installation folder.
1.0.1.4	30/01/2019	HD21ABE17 monitor bug fixing.
1.0.1.3	20/12/2018	Bugs fixing.
1.0.1.2	13/11/2018	Regenerated setup file.
1.0.1.1	09/11/2018	Added new functionalities to comply with ISO7243:2017 standard. Corrected a calculation inaccuracy of the ISO7933 PHS index.
0.1.6.7	17/10/2018	Added setting of CAV parameter for WBGT correction according to ISO 7243:2017 standard.